



EPIDURAL ANESTHESIA SIMULATOR

INSTRUCTIONS FOR CARE

1. Use ONLY needles with sharp bevel tip and diameter $\leq 22G$ to preserve the life-span of the ULTRASOUND SIMULATOR
2. Using needles of large caliber may create permanent damage, needle insertion track marks and require replacement of the tissue SIMULATOR insert sooner
3. **When practicing injection:**
 - a. Use only sterile water
 - b. Wipe off and dry the SIMULATOR before storing



Warranty

Your **ULTRASOUND SIMULATOR** includes 1-Year limited warranty against manufacturing and material defects. Visit www.NYSORASIMULATORS.com/support for further information to register your SIMULATOR and activate the warranty


Model: **MS2-EPI**

NYSORA 
SIMULATORS.COM

LEGEND:

SP - Spinous Process
AP - Articular Process

PL006 RevB



SIMULATOR IN USE

NYSORA  **SIMULATORS**TM

**PROXIMAL
EPIDURAL**